



NAILING SYSTEM

S.H.PITKAR
ORTHOTOOLS PVT. LTD.

Gat 36 & 37, RMK Industrial Park II, Talegaon MIDC Rd,
Jadhavwadi, Navalakh Umbare, Pune - 410507.
Mobile: +9175984730/31; Email: info@pitkar.net

ISO 13485 and Indian FDA approved company
CDSCO Approved
www.pitkar.com

Universal Locking Nail Titanium



Titanium Tibia Nail

Part No	Length	Dia
A 410008270	270 mm	8 mm
A 410008285	285 mm	8 mm
A 410008300	300 mm	8 mm
A 410008315	315 mm	8 mm
A 410008330	330 mm	8 mm
A 410008345	345 mm	8 mm
A 410008360	360 mm	8 mm
A 410008375	375 mm	8 mm
A 410008390	390 mm	8 mm
A 410009270	270 mm	9 mm
A 410009285	285 mm	9 mm
A 410009300	300 mm	9 mm
A 410009315	315 mm	9 mm
A 410009330	330 mm	9 mm
A 410009345	345 mm	9 mm
A 410009360	360 mm	9 mm
A 410009375	375 mm	9 mm
A 410009390	390 mm	9 mm
A 410010270	270 mm	10 mm
A 410010285	285 mm	10 mm
A 410010300	300 mm	10 mm
A 410010315	315 mm	10 mm
A 410010330	330 mm	10 mm
A 410010345	345 mm	10 mm
A 410010360	360 mm	10 mm
A 410010375	375 mm	10 mm
A 410010390	390 mm	10 mm
A 410011270	270 mm	11 mm
A 410011285	285 mm	11 mm
A 410011300	300 mm	11 mm
A 410011315	315 mm	11 mm
A 410011330	330 mm	11 mm
A 410011345	345 mm	11 mm
A 410011360	360 mm	11 mm
A 410011375	375 mm	11 mm
A 410011390	390 mm	11 mm

Pitkar

Pitkar

Pitkar

Titanium Femur Nail

Part No	Length	Dia
A 409009320	320 mm	9 mm
A 409009340	340 mm	9 mm
A 409009360	360 mm	9 mm
A 409009380	380 mm	9 mm
A 409009400	400 mm	9 mm
A 409009420	420 mm	9 mm
A 409009440	440 mm	9 mm
A 409009460	460 mm	9 mm
A 409009480	480 mm	9 mm
A 409010320	320 mm	10 mm
A 409010340	340 mm	10 mm
A 409010360	360 mm	10 mm
A 409010380	380 mm	10 mm
A 409010400	400 mm	10 mm
A 409010420	420 mm	10 mm
A 409010440	440 mm	10 mm
A 409010460	460 mm	10 mm
A 409010480	480 mm	10 mm
A 409011320	320 mm	11 mm
A 409011340	340 mm	11 mm
A 409011360	360 mm	11 mm
A 409011380	380 mm	11 mm
A 409011400	400 mm	11 mm
A 409011420	420 mm	11 mm
A 409011440	440 mm	11 mm
A 409011460	460 mm	11 mm
A 409011480	480 mm	11 mm
A 409012320	320 mm	12 mm
A 409012340	340 mm	12 mm
A 409012360	360 mm	12 mm
A 409012380	380 mm	12 mm
A 409012400	400 mm	12 mm
A 409012420	420 mm	12 mm
A 409012440	440 mm	12 mm
A 409012460	460 mm	12 mm
A 409012480	480 mm	12 mm

Pitkar

InterLocking Screw Titanium



Screw 4.9 mm

Part No	Length
A 408000030	30 mm
A 408000032	32 mm
A 408000034	34 mm
A 408000036	36 mm
A 408000038	38 mm
A 408000040	40 mm
A 408000042	42 mm
A 408000044	44 mm
A 408000046	46 mm
A 408000048	48 mm
A 408000050	50 mm
A 408000052	52 mm
A 408000054	54 mm
A 408000056	56 mm
A 408000058	58 mm
A 408000060	60 mm
A 408000065	65 mm
A 408000070	70 mm
A 408000075	75 mm
A 408000080	80 mm
A 408000085	85 mm
A 408000090	90 mm

Screw 3.9 mm

Part No	Length
A 411000030	30 mm
A 411000032	32 mm
A 411000034	34 mm
A 411000036	36 mm
A 411000038	38 mm
A 411000040	40 mm
A 411000042	42 mm
A 411000044	44 mm
A 411000046	46 mm

Pitkar

Pitkar

Pitkar

UNT Nail System



UNT Nail Tibia

Part No	Length	Dia
AI29708270	270 mm	8 mm
AI29708285	285 mm	8 mm
AI29708300	300 mm	8 mm
AI29708315	315 mm	8 mm
AI29708330	330 mm	8 mm
AI29708345	345 mm	8 mm
AI29708360	360 mm	8 mm
AI29708375	375 mm	8 mm
AI29708390	390 mm	8 mm
AI29709270	270 mm	9 mm
AI29709285	285 mm	9 mm
AI29709300	300 mm	9 mm
AI29709315	315 mm	9 mm
AI29709330	330 mm	9 mm
AI29709345	345 mm	9 mm
AI29709360	360 mm	9 mm
AI29709375	375 mm	9 mm
AI29709390	390 mm	9 mm
AI29710270	270 mm	10 mm
AI29710285	285 mm	10 mm
AI29710300	300 mm	10 mm
AI29710315	315 mm	10 mm
AI29710330	330 mm	10 mm
AI29710345	345 mm	10 mm
AI29710360	360 mm	10 mm
AI29710375	375 mm	10 mm
AI29710390	390 mm	10 mm
AI29711270	270 mm	11 mm
AI29711285	285 mm	11 mm
AI29711300	300 mm	11 mm
AI29711315	315 mm	11 mm
AI29711330	330 mm	11 mm
AI29711345	345 mm	11 mm
AI29711360	360 mm	11 mm
AI29711375	375 mm	11 mm
AI29711390	390 mm	11 mm

Pitkar

Pitkar



UNT Nail Femur

Part No	Length	Dia
AI29609320	320 mm	9 mm
AI29609340	340 mm	9 mm
AI29609360	360 mm	9 mm
AI29609380	380 mm	9 mm
AI29609400	400 mm	9 mm
AI29609420	420 mm	9 mm
AI29609440	440 mm	9 mm
AI29609460	460 mm	9 mm
AI29609480	480 mm	9 mm
AI29610320	320 mm	10 mm
AI29610340	340 mm	10 mm
AI29610360	360 mm	10 mm
AI29610380	380 mm	10 mm
AI29610400	400 mm	10 mm
AI29610420	420 mm	10 mm
AI29610440	440 mm	10 mm
AI29610460	460 mm	10 mm
AI29610480	480 mm	10 mm
AI29611320	320 mm	11 mm
AI29611340	340 mm	11 mm
AI29611360	360 mm	11 mm
AI29611380	380 mm	11 mm
AI29611400	400 mm	11 mm
AI29611420	420 mm	11 mm
AI29611440	440 mm	11 mm
AI29611460	460 mm	11 mm
AI29611480	480 mm	11 mm
AI29612320	320 mm	12 mm
AI29612340	340 mm	12 mm
AI29612360	360 mm	12 mm
AI29612380	380 mm	12 mm
AI29612400	400 mm	12 mm
AI29612420	420 mm	12 mm
AI29612440	440 mm	12 mm
AI29612460	460 mm	12 mm
AI29612480	480 mm	12 mm

Pitkar

Pitkar

Interlocking Screw Titanium



Screw 4.9 mm

Part No	Length
AI29800030	30 mm
AI29800032	32 mm
AI29800034	34 mm
AI29800036	36 mm
AI29800038	38 mm
AI29800040	40 mm
AI29800042	42 mm
AI29800044	44 mm
AI29800046	46 mm
AI29800048	48 mm
AI29800050	50 mm
AI29800052	52 mm
AI29800054	54 mm
AI29800056	56 mm
AI29800058	58 mm
AI29800060	60 mm
AI29800065	65 mm
AI29800070	70 mm
AI29800075	75 mm
AI29800080	80 mm
AI29800085	85 mm
AI29800090	90 mm

Screw 3.9 mm

Part No	Length
A 411000030	30 mm
A 411000032	32 mm
A 411000034	34 mm
A 411000036	36 mm
A 411000038	38 mm
A 411000040	40 mm
A 411000042	42 mm
A 411000044	44 mm
A 411000046	46 mm



Pitkar

Pitkar

Pitkar

- A 411000048 48 mm
- A 411000050 50 mm
- A 411000052 52 mm
- A 411000054 54 mm
- A 411000056 56 mm
- A 411000058 58 mm
- A 411000060 60 mm
- A 411000062 62 mm
- A 411000064 64 mm
- A 411000066 66 mm
- A 411000068 68 mm
- A 411000070 70 mm

Instrument Set for UNT Tibia Nail System
Part No. : AA 01110000

Guide Wire 2.5 mm
Part No. A 638000000



Depth Gauge
Part No. A 647001000



Extractor Hammer
Part No. A 651000000



Interlocking Awl
Part No. A 661001000



Combination Wrench 13 mm
Part No. D 502131000



Sleeve 11/8 mm
Part No. A 948000000



Drill Sleeve 8/4.3 mm
Part No. A 949000000



Drill Sleeve 8/3.3 mm
Part No. A 950000000



Trocar 8 mm
Part No. A 951000000



Screw Driver Hexagonal 3.5 A/F
Part No. A 952001000



Extractor Rod
Part No. A 953000000



Extractor Rod Handle
Part No. A 954000000



Drill Bit 3.2 mm
270 mm Extra Long
Part No. A 955027032



3.9 mm Screw - To be used for Distal Locking 8mm, 9mm of Universal Tibia & 9mm Universal Femur Nail.

4.9 mm Screw - To be used for
 1) Proximal Locking for Universal Tibia & Femur all diameter nails, Sol-T Nail and Yt Nail.
 2) Distal Locking of TFN all diameter nails.



Cap Screw for UNT Nail Tibia
Part No. AI29500000



Cap Screw for UNT Nail Femur
Part No. AI29400000



**Drill Bit 4.2 mm
270 mm Extra Long**
Part No.

A 955027042



Universal Locking Jig for Tibia
Part No.

AA 00701000



**Coupling Bolt for Universal
Locking Nail Jig**
Part No.

A 918000000

Impactor Attachment for Univerasl JIG
Part No.

AA07701000



Container for UNT Tibia Instrument Set
Part No.

H 935001000

**Instrument Set for UNT Distal
Tibia Nail System**

Part No : AA 01310000



Guide Wire 2.5 mm
Part No.

A 638000000



Depth Gauge
Part No.

A 647001000



Extractor Hammer
Part No.

A 651000000



Interlocking Awl
Part No.

A 661001000



Combination Wrench 13 mm
Part No.

D 502131000



Sleeve 11/8 mm
Part No.

A 948000000



Drill Sleeve 8/4.3 mm
Part No

A 949000000



Drill Sleeve 8/3.3 mm
Part No

A 950000000



Trocar 8 mm
Part No

A 951000000



Screw Driver Hexagonal 3.5 A/F
Part No.

A 952001000



Extractor Rod
Part No.

A 953000000



Extractor Rod Handle
Part No.

A 954000000



**Drill Bit 3.2 mm
270 mm Extra Long**
Part No.

A 955027032



**Drill Bit 4.2 mm
270 mm Extra Long**
Part No.

A 955027042



Proximal Distal Jig for Tibia
Part No. AA 00901000

Instrument Set for UNT Femur Nail System
Part No : AA 01210000



Coupling Bolt for UNT Nail Jig
Part No. A 918000000



Guide Wire 2.5 mm
Part No. A 638000000



Drill Sleeve 8/3.3 mm
Part No A 950000000



Container for UNT Distal Tibia Instrument Set
Part No. H 937001000



Depth Gauge
Part No. A 647001000



Trocar 8 mm
Part No A 951000000



Extractor Hammer
Part No. A 651000000



Screw Driver Hexagonal 3.5 A/F
Part No. A 952001000



Extractor Rod
Part No. A 953000000



Interlocking Awl
Part No. A 661001000



Extractor Rod
Part No. A 954000000



Combination Wrench 13 mm
Part No. D 502131000



Extractor Rod Handle
Part No. A 954000000



Sleeve 11/8 mm
Part No. A 948000000



Drill Sleeve 8/4.3 mm
Part No A 949000000



**Drill Bit 4.2 mm
270 mm Extra Long**
Part No.

A 955027042



Universal Locking Jig for Femur
Part No.

AA 00801000



Coupling Bolt for UNT Nail Jig
Part No.

A 918000000

**Impactor Attachment for
Universal Jig**
Part No.

AA07701000



**Container for UNT Femur Instrument
Set**
Part No.

H 936001000

**Instrument Set for UNT Distal
Femur Nail System**

Part No : AA 01410000



Guide Wire 2.5 mm
Part No.

A 638000000



Depth Gauge
Part No.

A 647001000



Extractor Hammer
Part No.

A 651000000



Interlocking Awl
Part No.

A 661001000



Combination Wrench 13 mm
Part No.

D 502131000



Sleeve 1 1/8 mm
Part No.

A 948000000



Drill Sleeve 8/4.3 mm
Part No

A 949000000



Drill Sleeve 8/3.3 mm
Part No

A 950000000



Trocar 8 mm
Part No

A 951000000



Screw Driver Hexagonal 3.5 A/F
Part No.

A 952001000



Extractor Rod
Part No.

A 953000000



Extractor Rod Handle
Part No.

A 954000000



**Drill Bit 3.2 mm
270 mm Extra Long**
Part No.

A 955027032



**Drill Bit 4.2 mm
270 mm Extra Long**
Part No.

A 955027042



Proximal Distal Jig for Femur
Part No.

AA 01001000



**Coupling Bolt for Universal
Locking Nail Jig**
Part No.

A 918000000



**Container For UNT Distal Femur
Instrument Set**
Part No.

H 938001000

SOL - T Nail System Titanium



SOL-T- Tibia Nail System

Part No	Length	Dia
A 447008270	270 mm	8mm
A 447008285	285 mm	8mm
A 447008300	300 mm	8 mm
A 447008315	315 mm	8 mm
A 447008330	330 mm	8 mm
A 447008345	345 mm	8 mm
A 447008360	360 mm	8 mm
A 447008375	375 mm	8 mm
A 447009270	270 mm	9 mm
A 447009285	285 mm	9 mm
A 447009300	300 mm	9 mm
A 447009315	315 mm	9 mm
A 447009330	330 mm	9 mm
A 447009345	345 mm	9 mm
A 447009360	360 mm	9 mm
A 447009375	375 mm	9 mm
A 447010270	270 mm	10 mm
A 447010285	285 mm	10 mm
A 447010300	300 mm	10 mm
A 447010315	315 mm	10 mm
A 447010330	330 mm	10 mm
A 447010345	345 mm	10 mm
A 447010360	360 mm	10 mm
A 447010375	375 mm	10 mm

SOL-T Femur Nail System

Part No	L	Right	Dia (mm)	
Left	Right			
A 445009320	A 446009320	320	9	
A 445009340	A 446009340	340	9	
A 445009360	A 446009360	360	9	
A 445009380	A 446009380	380	9	
A 445009400	A 446009400	400	9	
A 445009420	A 446009420	420	9	
A 445009440	A 446009440	440	9	
A 445009460	A 446009460	460	9	
A 445009480	A 446009480	480	9	
A 445010320	A 446010320	320	10	
A 445010340	A 446010340	340	10	
A 445010360	A 446010360	360	10	
A 445010380	A 446010380	380	10	
A 445010400	A 446010400	400	10	
A 445010420	A 446010420	420	10	
A 445010440	A 446010440	440	10	
A 445010460	A 446010460	460	10	
A 445010480	A 446010480	480	10	
A 445011320	A 446011320	320	11	
A 445011340	A 446011340	340	11	
A 445011360	A 446011360	360	11	
A 445011380	A 446011380	380	11	
A 445011400	A 446011400	400	11	
A 445011420	A 446011420	420	11	
A 445011440	A 446011440	440	11	
A 445011460	A 446011460	460	11	
A 445011480	A 446011480	480	11	



Cervical Screws Dia 6.5 mm

Part No

- A 448000060
- A 448000065
- A 448000070
- A 448000075
- A 448000080
- A 448000085
- A 448000090
- A 448000095
- A 448000100
- A 448000105
- A 448000110

Length

- 60 mm
- 65 mm
- 70 mm
- 75 mm
- 80 mm
- 85 mm
- 90 mm
- 95 mm
- 100 mm
- 105 mm
- 110 mm

SOL-T Tibia Instruments Set

Part No. A 749000000



Jig for straight hole

Part No. A 675001000



Jig for 45deg. hole

Part No. A 676001000



Tissue Protection Sleeve 10.5 mm

Part No. A 816000000



Coupling bolt

Part No. A 677000000



Reamer 10.5 mm

Part No. A 817000000



2.5*1000 mm Unbided Guide Wire

Part No. A 638000000



Medullary Tube Polyamide

Part No. A 653001000



K-awl

Part No. A 661000000



8 mm Trocar

Part No. A 951000000



Tissue Protector

Part No. A 654000000



Depth Gauge

Part No. A 647001000



Screw Driver Hexagonal 3.5 mm A/F

Part No. A 952001000



Extractor rod

Part No. A 953000000



Extractor Hammer

Part No. A 651000000



Extractor Handle

Part No. A 954000000



Combination Wrench 13 mm

Part No. D 502131000



T Handle

Part No. A 945001000



Protection Sleeve 11/8mm

Part No. A 948000000



8/3.5*137 mm Drill Guide

Part No. A 963001000



8/3*137 mm Drill Guide

Part No. A 871001000



Drill Bit 3.5 x 250 mm

Part No. A 964025035



Drill Bit 3 x 250 mm

Part No. A 866000000



Universal Spanner A/F 13mm

Part No. A 965001000



Containers for SOL-T System

Part No.

Implants H 846001000

Instruments H 812101000

SOL-T Femur Instruments Set

Part No. A 750000000



Jig for straight hole

Part No. A 706001000



Measuring Device

Part No. A 699000000



Tissue Protection Sleeve 13.5 mm

Part No. A 815001000



Coupling bolt

Part No. A 918000000



13.5 mm Quick Change Reamer (300mm)

Part No. A 814001000



2.5*1000 mm Unbided Guide Wire

Part No. A 638000000



Medullary Tube Polyamide

Part No. A 653001000



K-awl

Part No. A 661000000



8 mm Trocar

Part No. A 951000000



Tissue Protector

Part No. A 654000000



Depth Gauge

Part No. A 647001000



Screw Driver Hexagonal 3.5 mm A/F

Part No. A 952001000



Extractor rod

Part No. A 953000000



Extractor Hammer

Part No. A 651000000



Extractor Handle

Part No. A 954000000



Combination Wrench 13 mm

Part No. D 502131000



T Handle

Part No. A 945001000



Protection Sleeve 11/8mm

Part No. A 948000000



8/3.5*137 mm Drill Guide

Part No. A 963001000



8/3*137 mm Drill Guide

Part No. A 871001000



Drill Bit 3.5 x 250 mm

Part No. A 964025035



Drill Bit 3 x 250 mm

Part No. A 866000000



Universal Spanner A/F 13mm

Part No. A 965001000



Protection Sleeve 10/8 X 160 mm

Part No. A 962000000



Sleeve 8/2 x 160mm
Part No. A 704001000



6.5mm step reamer
Part No. A 701001000



**3.5mm Hex Screw Driver
Cannulated**
Part No. A 700001000



Guide Wire 2x390mm
Part No. A 708000000



**Supine drive for
TFN nail**
Part No. AA05700000



**Bolt locating attachment
for TFN nail**
Part No. AA05901000



**Hex nut for bolt locating
attachment**
Part No. AA05800000



Pilot Reamer Ø 6mm
Part No. AA07200000



Containers for SOL-T System
Part No.
Implants H 846001000
Instruments H 812201000

SOL-T Proximal Femoral Device



SOL-T PFD Nail Cannulated
Part No
A 443009242 Dia 9 mm
A 443010242 10 mm
A 443011242 11 mm



SOL-T PFD Solid Nail
Part No. Dia
A 456009130 9 mm
A 456010130 10 mm
A 456011130 11 mm



SOL-T PFD Long Nail System

Part No	L	Right	Dia	mm
A 365009340	A 366009340	340	9	
A 365009380	A 366009380	380	9	
A 365009400	A 366009400	400	9	
A 365009420	A 366009420	420	9	
A 365010340	A 366010340	340	10	
A 365010380	A 366010380	380	10	
A 365010400	A 366010400	400	10	
A 365010420	A 366010420	420	10	
A 365011340	A 366011340	340	11	
A 365011380	A 366011380	380	11	
A 365011400	A 366011400	400	11	
A 365011420	A 366011420	420	11	



Hip Pin 6 mm

Part No	Length
A 43000055	55 mm
A 43000060	60 mm
A 43000065	65 mm
A 43000070	70 mm
A 43000075	75 mm
A 43000080	80 mm
A 43000085	85 mm
A 43000090	90 mm
A 43000095	95 mm
A 43000100	100 mm



Hip Bolt 9 mm

Part No	Length
A 43100070	70 mm
A 43100075	75 mm
A 43100080	80 mm
A 43100085	85 mm
A 43100090	90 mm
A 43100095	95 mm
A 43100100	100 mm
A 43100105	105 mm
A 43100110	110 mm
A 43100115	115 mm

Cap Screw for SOL-T PFD Nail

Part No. A358900000

**SOL-T Proximal Femoral Device
Instruments Set- A 751000000**

Ram

Part No. A 696000000



Ram Sleeve

Part No. A 697000000



Ram Rod

Part No. A 698001000



Combination Wrench 11 mm

Part No. A 702000000



Depth Gauge

Part No. A 709001000



**Universal Spanner for
Tibia/Femur Sol-T Set**

Part No. A 818001000



Sleeve 12 x 10 mm

Part No. A 725001000



Sleeve 10 x 2.5 mm

Part No. A 726001000



Trocar 2.5 mm for 12x10 mm Sleeve

Part No. A 733001000



Sleeve 9 x 7 mm

Part No. A 727001000



Sleeve 7 x 2.5 mm

Part No. A 728001000



**Trocar 2.5 mm for
9 x 7 sleeve**

Part No. A 748001000



Sleeve 18 x 15 mm

Part No. A 729001000





Sleeve 15 x 2.5 mm
Part No. A 730001000



**Cannulated Screw Driver
A/F 4x400mm**
Part No. A 736001000



Measuring Device
Part No. A 745000000



Trocar 2.5 mm for 18 x 15 Sleeve
Part No. A 743000000



Sol-T PFD Jig
Part No. A 840001000



Guide Wire 2.5x440 mm
Part No. A 744000000



Tissue Protection Sleeve 10/8 mm
Part No. A 685001000



**T Quick Coupling
A/F 12mm Drill**
Part No. A 747001000



Drill Guide 8x3.5mm
Part No. A 872001000



6mm Step Reamer
Part No. A 738000000



Drill Bit 3.5 x 350 mm
Part No. A 865000000



9mm Step Reamer
Part No. A 739000000



Guide Wire 2.5x 350 mm
Part No. A 869000000



Trocar 3.5 mm
Part No. A 874001000



Reamer 15mm
Part No. A 741000000



Supine drive for TFN nail
Part No. AA05700000



Wrench for Hip bolt
Part No. A 735001000



Coupling Bolt
Part No. A 742000000



Screw Driver Solid 3.5 mm A/F
Part No. A 746001000



**Bolt locating attachment for
TFN nail**
Part No. AA05901000



Hex nut for bolt locating attachment A/F 11mm

Part No. AA07700000



Pilot Reamer Ø 6mm

Part No. AA07800000



Container for Sol-T Proximal Femoral Device

Description	Part No.
Implants	H 844001000
Instruments	H 820001000

Supracondylar Nail



Part No	Length	Dia
A 367010150	150 mm	10 mm
A 367010250	250 mm	10 mm
A 367010300	300 mm	10 mm
A 367010350	350 mm	10 mm
A 367011150	150 mm	11 mm
A 367011250	250 mm	11 mm
A 367011300	300 mm	11 mm
A 367011350	350 mm	11 mm
A 367012150	150 mm	12 mm
A 367012250	250 mm	12 mm
A 367012300	300 mm	12 mm
A 367012350	350 mm	12 mm

A 432000025	25 mm
A 432000030	30 mm
A 432000035	35 mm
A 432000040	40 mm
A 432000045	45 mm
A 432000050	50 mm
A 432000055	55 mm
A 432000060	60 mm
A 432000065	65 mm
A 432000070	70 mm
A 432000075	75 mm
A 432000080	80 mm
A 432000085	85 mm
A 432000090	90 mm
A 432000095	95 mm
A 432000100	100 mm

Interlocking Screw



Interlocking Screw 5 mm

Part No	Length
A 433000035	35 mm
A 433000040	40 mm
A 433000045	45 mm
A 433000050	50 mm
A 433000055	55 mm
A 433000060	60 mm
A 433000065	65 mm
A 433000070	70 mm
A 433000075	75 mm

Interlocking Screw 3.8 mm

Part No	Length
A 449000020	20 mm
A 449000025	25 mm
A 449000030	30 mm
A 449000035	35 mm
A 449000040	40 mm
A 449000045	45 mm
A 449000050	50 mm
A 449000055	55 mm
A 449000060	60 mm
A 449000065	65 mm
A 449000070	70 mm
A 449000075	75 mm



Interlocking Screw 4.5 mm

Part No	Length
A 432000020	20 mm

Trocantric Femoral Nail



Interlocking Trocantric Nail

Part No.	Length	Dia
A 450130009	130	9 mm
A 450130010	130	10 mm
A 450130011	130	11 mm
A 450135009	135	9 mm
A 450135010	135	10 mm
A 450135011	135	11 mm



Anti Rotational Screw 6.4 mm

Part No.	Dia
A 451064065	65 mm
A 451064070	70 mm
A 451064075	75 mm
A 451064080	80 mm
A 451064085	85 mm
A 451064090	90 mm
A 451064095	95 mm
A 451064100	100 mm
A 451064105	105 mm
A 451064110	110 mm



Lag Screw 8 mm

Part No.	Dia
A 452008065	65 mm
A 452008070	70 mm
A 452008075	75 mm
A 452008080	80 mm
A 452008085	85 mm
A 452008090	90 mm
A 452008095	95 mm
A 452008100	100 mm
A 452008105	105 mm

Trocantric Femoral Nail Instrument Set

Part No. A 921010000



8mm Step Reamer for TFN Cannulated

Part No. A 922000000



6.4mm Step Reamer for TFN Cannulated

Part No. A 923000000



Drill Bit

Dia 4X250 mm Part No. A 924000000



Jig

Part No. A 926001000
A 925000100
130°
135°



Bolt for Trocantric Nail

Part No. A 927000000

Protection Guide Sleeve
12 X 10mm

Part No. A 928000000

Protection Guide Sleeve
10 X 8 mm

Part No. A 930000000



Drill Sleeve 8 X 4mm Long

Part No. A 931000000



Guide Wire Sleeve

Part No. A 933000000
A 932000000
8X2.5 mm
10X2.5 mm



TFN Trocar Sleeve 4mm

Part No. A 934001000

Hexagonal Screw Driver
6.4mm Cannulated

Part No. A 938001000



**Hex Screw Driver
4.5 mm Extra Long
Part No.**

A 937001000



**Combination Wrench
Across Flat 10 mm
Part No.**

A 630000000



**Extractor Rod Handle
Part No.**

A 954000000



**T Handle
Part No.**

A 945001000



**Extractor Hammer
Part No.**

A 651000000



**Measuring Device
Part No.**

A 994000000



**Calibrated Guide Wire
Threaded 2.5 X 390 mm
Part No.**

A 947000000



**Extractor Rod
Part No.**

A 953000000



**Container
Part No.**

A 913001000



Titanium Yt Nail 180mm

Part No.	Dia (mm)	Deg
A 479009180*	9	120°
A 478009180*	9	125°
A 475009180*	9	130°
A 479010180	10	120°
A 478010180	10	125°
A 475010180	10	130°
A 479011180	11	120°
A 478011180	11	125°
A 475011180	11	130°

Titanium Yt Nail 210mm

Part No.	Dia (mm)	Deg
A 479009210*	9	120°
A 478009210*	9	125°
A 475009210*	9	130°
A 479010210*	10	120°
A 478010210*	10	125°
A 475010210*	10	130°
A 479011210	11	120°
A 478011210	11	125°
A 475011210	11	130°

Titanium Yt Nail 240mm

Part No.	Dia (mm)	Deg
A 479009240*	9	120°
A 478009240*	9	125°
A 475009240*	9	130°
A 479010240*	10	120°
A 478010240*	10	125°
A 475010240*	10	130°

A 479011240*	11	120°
A 478011240*	11	125°
A 475011240*	11	130°

Yt Long Nail System

Yt Long Nail 360mm Dia 10mm

Part No.	Left	Right	Deg
A 484010360	A 483010360		120°
A 486010360	A 485010360		125°
A 488010360	A 487010360		130°

Yt Long Nail 380mm Dia 10mm

Part No.	Left	Right	Deg
A 484010380	A 483010380		120°
A 486010380	A 485010380		125°
A 488010380	A 487010380		130°

Yt Long Nail 400mm Dia 10mm

Part No.	Left	Right	Deg
A 484010400	A 483010400		120°
A 486010400	A 485010400		125°
A 488010400	A 487010400		130°

Yt Long Nail 420mm Dia 10mm

Part No.	Left	Right	Deg
A 484010420	A 483010420		120°
A 486010420	A 485010420		125°
A 488010420	A 487010420		130°

* Available only on orders



Yt nail Screw

Part No.	Dia	Length
A 473010070	10	70
A 473010075	10	75
A 473010080	10	80
A 473010085	10	85
A 473010090	10	90
A 473010095	10	95
A 473010100	10	100
A 473010105	10	105
A 473010110	10	110
A 473010115	10	115
A 473010120	10	120



Interlocking Screw 4.9mm

Part No.	Length
AI 29800024	24
AI 29800026	26
AI 29800028	28
AI 29800030	30
AI 29800032	32
AI 29800034	34
AI 29800036	36
AI 29800038	38
AI 29800040	40
AI 29800042	42
AI 29800044	44
AI 29800046	46
AI 29800048	48
AI 29800050	50
AI 29800052	52
AI 29800054	54
AI 29800056	56
AI 29800058	58
AI 29800060	60
AI 29800065	65
AI 29800070	70
AI 29800075	75
AI 29800080	80
AI 29800085	85
AI 29800090	90

Instruments for Yt Nail



Set Screw for Yt Nail
Part No. A 474001000



Cap Screw for Yt Nail
Part No. A 477000000

Yt nail Instrument Set
Part No. AA 04010000



Measuring Device
Part No. AA 02400000



Protection Sleeve 11 x 8 mm
Part No. AA 02500000



Drill Sleeve 8x4.2mm
Part No. AA 02600000



Trocar 8mm
Part No. AA 02700000



Reamer 15mm
Part No. AA 02800000



Reamer 10mm
Part No. AA 02900000



Yt nail Jig
Part No. AA 03001000



Jig Cover for 120°
Part No. AA 03001200



Jig Cover for 125°
Part No. AA 03001300



Jig Cover for 130°
Part No. AA 03001400



Coupling Bolt
Part No. AA 03100000



Universal Spanner
Part No. AA 03201000



Introducer for Yt nail Screw
Part No. AA 03301000



Guide Wire 2.5x382mm
Part No. AA 03400000



Guide Wire 2.5x500mm
Part No. AA 03500000



Protection Sleeve 14 x 10mm
Part No. AA 03600000



Guide Wire Sleeve 10x2.5mm
Part No. AA 03700000



Impactor

Part No. AA 03801000



Sleeve 18x15mm

Part No. A 729001000



Sleeve 15x2.5mm

Part No. A 730001000



Trocar 2.5mm

Part No. A 743001000



Screw Driver for Yt nail

Part No. AA 04101000



Extra Long Drill Bit 4.2 x 270 mm

Part No. A 955027042



T Handle

Part No. A 945001000



Tommy Bar for Yt nail

Part No. AA 05000000



Hex Extension Sleeve

Part No. AA 05300000




Depth gauge

Part No. A647001000



Cannulated Curved Awl

Part No. AA 05401000



Container for Yt nail Instrument Set

Part No. H 943001000

PFN Next Gen Nail System



PFN Next Gen Nail 180mm

Part No.	Dia (mm)	Deg
AI28809180	9	130°
AI28810180	10	130°
AI28811180	11	130°
AI28812180	12	130°
AI28813180*	13	130°
AI28909180	9	125°
AI28910180	10	125°
AI28911180	11	125°
AI28912180	12	125°
AI28913180*	13	125°



PFN Next Gen Nail 210mm

Part No.	Dia	Deg
AI28809210	9	130°
AI28810210	10	130°
AI28811210	11	130°
AI28812210	12	130°
AI28813210*	13	130°
AI28909210	9	125°
AI28910210	10	125°
AI28911210	11	125°
AI28912210	12	125°
AI28913210*	13	125°

* Available only on orders



PFN Next Gen Nail 240mm

Part No.	Dia (mm)	Deg
AI28809240	9	130°
AI28810240	10	130°
AI28811240	11	130°
AI28812240	12	130°
AI28813240*	13	130°
AI28909240	9	125°
AI28910240	10	125°
AI28911240	11	125°
AI28912240	12	125°
AI28913240*	13	125°



Long PFN Next Gen Nail Ø 9 X 130°

Part No. Left	Part No. Right	Dia (mm)
AI29009320	AI29109320	9
AI29009340	AI29109340	9
AI29009360*	AI29109360*	9
AI29009380*	AI29109380*	9
AI29009400*	AI29109400*	9
AI29009420*	AI29109420*	9

Long PFN Next Gen Nail Ø 10 X 130°

Part No. Left	Part No. Right	Dia (mm)
AI29010320	AI29110320	10
AI29010340	AI29110340	10
AI29010360*	AI29110360*	10
AI29010380	AI29110380	10
AI29010400*	AI29110400*	10
AI29010420	AI29110420	10

Long PFN Next Gen Nail Ø 11 X 130°

Part No. Left	Part No. Right	Dia (mm)
AI29011320*	AI29111320*	11
AI29011340*	AI29111340*	11
AI29011360*	AI29111360*	11
AI29011380*	AI29111380*	11
AI29011400*	AI29111400*	11
AI29011420*	AI29111420*	11

Long PFN Next Gen Nail Ø 12 X 130°

Part No. Left	Part No. Right	Dia (mm)
AI29012320*	AI29112320*	12
AI29012340*	AI29112340*	12
AI29012360*	AI29112360*	12
AI29012380*	AI29112380*	12
AI29012400*	AI29112400*	12
AI29012420*	AI29112420*	12

Long PFN Next Gen Nail Ø 13 X 130°

Part No. Left	Part No. Right	Dia (mm)
AI29013320*	AI29113320*	13
AI29013340*	AI29113340*	13
AI29013360*	AI29113360*	13
AI29013380*	AI29113380*	13
AI29013400*	AI29113400*	13
AI29013420*	AI29113420*	13

Long PFN Next Gen Nail Ø 9 X 125°

Part No. Left	Part No. Right	Dia (mm)
AI29209320	AI29309320	9
AI29209340	AI29309340	9
AI29209360*	AI29309360*	9
AI29209380*	AI29309380*	9
AI29209400*	AI29309400*	9
AI29209420*	AI29309420*	9

Long PFN Next Gen Nail \varnothing 10 X 125°

Part No. Left	Part No. Right	Dia (mm)
AI29210320	AI29310320	10
AI29210340	AI29310340	10
AI29210360*	AI29310360*	10
AI29210380	AI29310380	10
AI29210400*	AI29310400*	10
AI29210420	AI29310420	10

Long PFN Next Gen Nail \varnothing 11 X 125°

Part No. Left	Part No. Right	Dia (mm)
AI29211320*	AI29311320*	11
AI29211340*	AI29311340*	11
AI29211360*	AI29311360*	11
AI29211380*	AI29311380*	11
AI29211400*	AI29311400*	11
AI29211420*	AI29311420*	11

Long PFN Next Gen Nail \varnothing 12 X 125°

Part No. Left	Part No. Right	Dia (mm)
AI29212320*	AI29312320*	12
AI29212340*	AI29312340*	12
AI29212360*	AI29312360*	12
AI29212380*	AI29312380*	12
AI29212400*	AI29312400*	12
AI29212420*	AI29312420*	12

Long PFN Next Gen Nail \varnothing 13 X 125°

Part No. Left	Part No. Right	Dia (mm)
AI29213320*	AI29313320*	13
AI29213340*	AI29313340*	13
AI29213360*	AI29313360*	13
AI29213380*	AI29313380*	13
AI29213400*	AI29313400*	13
AI29213420*	AI29313420*	13



*** SS Long PFNA-II Nail \varnothing 9 X 130°**

Part No. Left	Part No. Right	Dia (mm)
AI30509320	AI30609320	9
AI30509340	AI30609340	9
AI30509360	AI30609360	9
AI30509380	AI30609380	9
AI30509400	AI30609400	9
AI30509420	AI30609420	9

*** SS Long PFNA-II Nail \varnothing 10 X 130°**

Part No. Left	Part No. Right	Dia (mm)
AI30510320	AI30610320	10
AI30510340	AI30610340	10
AI30510360	AI30610360	10
AI30510380	AI30610380	10
AI30510400	AI30610400	10
AI30510420	AI30610420	10

*** SS Long PFNA-II Nail \varnothing 11 X 130°**

Part No. Left	Part No. Right	Dia (mm)
AI30511320	AI30611320	11
AI30511340	AI30611340	11
AI30511360	AI30611360	11
AI30511380	AI30611380	11
AI30511400	AI30611400	11
AI30511420	AI30611420	11

*** SS Long PFNA-II Nail \varnothing 12 X 130°**

Part No. Left	Part No. Right	Dia (mm)
AI30512320	AI30612320	12
AI30512340	AI30612340	12
AI30512360	AI30612360	12
AI30512380	AI30612380	12
AI30512400	AI30612400	12
AI30512420	AI30612420	12

*** SS Long PFNA-II Nail \varnothing 13 X 130°**

Part No. Left	Part No. Right	Dia (mm)
AI30513320	AI30613320	13
AI30513340	AI30613340	13
AI30513360	AI30613360	13
AI30513380	AI30613380	13
AI30513300	AI30613400	13
AI30513320	AI30613420	13

*** SS Long PFNA-II Nail \varnothing 9 X 125°**

Part No. Left	Part No. Right	Dia (mm)
AI30709320	AI30809320	9
AI30709340	AI30809340	9
AI30709360	AI30809360	9
AI30709380	AI30809380	9
AI30709400	AI30809400	9
AI30709420	AI30809420	9

*** SS Long PFNA-II Nail \varnothing 10 X 125°**

Part No. Left	Part No. Right	Dia (mm)
AI30710320	AI30810320	10
AI30710340	AI30810340	10
AI30710360	AI30810360	10
AI30710380	AI30810380	10
AI30710400	AI30810400	10
AI30710420	AI30810420	10

*** SS Long PFNA-II Nail \varnothing 11 X 125°**

Part No. Left	Part No. Right	Dia (mm)
AI30711320	AI30811320	11
AI30711340	AI30811340	11
AI30711360	AI30811360	11
AI30711380	AI30811380	11
AI30711400	AI30811400	11
AI30711420	AI30811420	11

* Available only on orders

*** SS Long PFNA-II Nail Ø 12 X 125°**

Part No.	Part No.	Dia
Left	Right	(mm)
AI30712320	AI30812320	12
AI30712340	AI30812340	12
AI30712360	AI30812360	12
AI30712380	AI30812380	12
AI30712400	AI30812400	12
AI30712420	AI30812420	12

*** SS Long PFNA-II Nail Ø 13 X 125°**

Part No.	Part No.	Dia
Left	Right	(mm)
AI30713320	AI30813320	13
AI30713340	AI30813340	13
AI30713360	AI30813360	13
AI30713380	AI30813380	13
AI30713400	AI30813400	13
AI30713420	AI30813420	13



PFN Next Gen Blade

Part No.	Dia	Length
	(mm)	(mm)
AI28701065*	10.5	65
AI28701070*	10.5	70
AI28701075	10.5	75
AI28701080	10.5	80
AI28701085	10.5	85
AI28701090	10.5	90
AI28701095	10.5	95
AI28701100	10.5	100
AI28701105	10.5	105
AI28701110	10.5	110
AI28701115	10.5	115
AI28701120	10.5	120



Interlocking Screw for PFN Next Gen Nail 4.9 mm Ti

Part No.	Length
A408000030	30
A408000032	32
A408000034	34
A408000036	36
A408000038	38
A408000040	40
A408000042	42
A408000044	44
A408000046	46
A408000048	48
A408000050	50
A408000052	52
A408000054	54
A408000056	56
A408000058	58
A408000060	60
A408000065	65
A408000070	70
A408000075	75
A408000080	80
A408000085	85
A408000090	90



Long Cap Screw for PFN Next Gen Nail

Part No.	
AI31300125	125°
AI31300130	130°



Cap Screw for PFN Next Gen Nail

Part No.	
AI30000000	

PFN Next Gen Nail Instrument Set



Cannulated Curved Awl

Part No.	
AA 05401000	



Sleeve 20 x 16.5mm

Part No.	
AA 06601000	



Sleeve 16.5 x 2.5mm

Part No.	
AA 06501000	



Trocar 2.5mm

Part No.	
AA 06401000	

* Available only on orders

Titanium Elastic Nail System



Elastic Nail

Part No.	Dia mm	Length mm
A 454020440	2.0	440
A 454025440	2.5	440
A 454030440	3.0	440
A 454035440	3.5	440
A 454040440	4.0	440

Instruments



Awl

Part No. A 996001000



Detachable Slide Hammer

Part No. A 997001000



Locking Plier

Part No. A 998001000



F Tool

Part No. A 999001000



Beveled Tamp

Part No. AA 00100000



Standard Tamp

Part No. AA 00300000



Spanner Wrench

Part No. AA 00200000



Inserter with Universal Chuck

Part No. AA 00401000

Container for
Elastic Nail
Instruments and
Implants Set

H 925001000

****Also Available On Request**



Double Drill Sleeve 4.5X3.2 mm

Part No. A 760001000



Drill Bit 3.2X125 mm

Part No. A 616012532



Drill Bit 4.5X125 mm

Part No. A 616012545



Drill Bit 2.8X150 mm

Part No. A 875000000

Guide Wire 2.5 x 500mm
Part No. AA03500000



Reamer 16.5mm
Part No. AA06300000



PFN Next Gen Jig
Part No. AA07101000



Jig Cover for 125°
Part No. AA07101300



Jig Cover for 130°
Part No. AA07101400



Coupling Bolt
Part No. AA03100000



Impactor
Part No. AA03801000



Tommy Bar for PFN Next Gen Nail
Part No. AA05000000



Universal Spanner
Part No. AA03201000



Hex Extension Sleeve
Part No. AA05300000



Protection Sleeve 15 x 11mm
Part No. AA06701000



Guide Wire Sleeve 11 x 3.2mm
Part No. AA06800000



Guide Wire 3.2 x 382mm
Part No. AA06100000



Measuring Device
Part No. AA07500000



Pilot Reamer 11mm
Part No. AA06900000



Step Reamer 11mm
Part No. AA07000000



T Handle
Part No. A945001000



PFN Next Gen Blade Introducer
Part No. AA06001000



Compression Instrument for PFN Next Gen Blade
Part No. AA07601000



Protection Sleeve 11 x 8mm
Part No. AA02500000



Trocar 8mm
Part No. AA02700000



Drill Sleeve 8 x 4.2mm
Part No. AA02600000



Extra Long Drill Bit 4.2 x 270mm
Part No. A955027042



Depth Gauge
Part No. A647001000



Screw Driver
Part No. AA04101000

Container for PFN Next Gen Instrument Set
Part No. H962001000

The Ultra Lock II Nails

Ultra Lock II Tibia Nail

- For Tibial mid shaft fractures both open and closed.
- For Pseudoarthrosis.
- Proximal 15 degree Bend-more anatomical.
- Dynamic Mediolateral Proximal Locking.
- Anterio Posterior Distal Locking
- Economical option of Nail.



Part No	Length (mm)	Dia (mm)
A 352008270	270 mm	8 mm
A 352008285	285 mm	8 mm
A 352008300	300 mm	8 mm
A 352008315	315 mm	8 mm
A 352008330	330 mm	8 mm
A 352008345	345 mm	8 mm
A 352008360	360 mm	8 mm
A 352008375	375 mm	8 mm
A 352008390	390 mm	8 mm
A 352009270	270 mm	9 mm
A 352009285	285 mm	9 mm
A 352009300	300 mm	9 mm
A 352009315	315 mm	9 mm
A 352009330	330 mm	9 mm
A 352009345	345 mm	9 mm
A 352009360	360 mm	9 mm
A 352009375	375 mm	9 mm
A 352009390	390 mm	9 mm
A 352010270	270 mm	10 mm
A 352010285	285 mm	10 mm
A 352010300	300 mm	10 mm
A 352010315	315 mm	10 mm
A 352010330	330 mm	10 mm
A 352010345	345 mm	10 mm
A 352010360	360 mm	10 mm
A 352010375	375 mm	10 mm
A 352010390	390 mm	10 mm

A 352011270	270 mm	11 mm
A 352011285	285 mm	11 mm
A 352011300	300 mm	11 mm
A 352011315	315 mm	11 mm
A 352011330	330 mm	11 mm
A 352011345	345 mm	11 mm
A 352011360	360 mm	11 mm
A 352011375	375 mm	11 mm
A 352011390	390 mm	11 mm

Ultra Lock II Femur Nail

- All mid shaft open & closed femoral fractures.
- Delayed unions.
- Universal (R or L), Anatomic Curve.
- Proximal Locking (greater trochanter to lesser trochanter) for improved thread purchase proximally.
- Transverse Distal Locking.
- Economical option of Nail.



Part No	Length	Dia
A 341009320	320 mm	9 mm
A 341009340	340 mm	9 mm
A 341009360	360 mm	9 mm
A 341009380	380 mm	9 mm
A 341009400	400 mm	9 mm
A 341009420	420 mm	9 mm
A 341009440	440 mm	9 mm
A 341009460	460 mm	9 mm
A 341009480	480 mm	9 mm
A 341010320	320 mm	10 mm
A 341010340	340 mm	10 mm
A 341010360	360 mm	10 mm
A 341010380	380 mm	10 mm
A 341010400	400 mm	10 mm
A 341010420	420 mm	10 mm
A 341010440	440 mm	10 mm
A 341010460	460 mm	10 mm

A 341010480	480 mm	10 mm
A 341011320	320 mm	11 mm
A 341011340	340 mm	11 mm
A 341011360	360 mm	11 mm
A 341011380	380 mm	11 mm
A 341011400	400 mm	11 mm
A 341011420	420 mm	11 mm
A 341011440	440 mm	11 mm
A 341011460	460 mm	11 mm
A 341011480	480 mm	11 mm
A 341012320	320 mm	12 mm
A 341012340	340 mm	12 mm
A 341012360	360 mm	12 mm
A 341012380	380 mm	12 mm
A 341012400	400 mm	12 mm
A 341012420	420 mm	12 mm
A 341012440	440 mm	12 mm
A 341012460	460 mm	12 mm
A 341012480	480 mm	12 mm

Ultra Lock II-Supracondylar Nail

- All Femoral supracondylar fractures
- Pathologic fractures.
- Failed distal femoral osteosynthesis and distal fractures.
- The 8 deg anterior bend makes easy insertion in the medullary canal.
- The nails are designed to permit the distal driving end to be counter sunk below the level of articular surface of femur.
- The first tow distal holes are 15 and 30 mm from the driving end.



Part No	Length	Dia
A 354010150	150 mm	10 mm
A 354010200	200 mm	10 mm
A 354010250	250 mm	10 mm
A 354011150	150 mm	11 mm
A 354011200	200 mm	11 mm
A 354011250	250 mm	11 mm

A 354012150	150 mm	12 mm
A 354012200	200 mm	12 mm
A 354012250	250 mm	12 mm
A 354013150	150 mm	13 mm
A 354013200	200 mm	13 mm
A 354013250	250 mm	13 mm

Instruments UltraLock II Nails

Instruments Set for Tibia Nails	AA 01710000
Instruments Set for Femur Nails	AA 01810000
Instruments Set for Distal Tibia Nails	AA 01910000
Instruments Set for Distal Femur Nails	AA 02010000
Instruments Set for Supra Nails	AA 02110000

SOL - T Nail System

SOL-T- Tibia Nail System

- The runner-shaped tip allows the implant to slide freely along the dorsal inner cortex.
- The special cross-section shape allows room for revascularisation and this favours the healing of the fracture.
- Can be used effectively for upper & lower one third tibial fractures & applications.



Part No	Length	Dia
A 350008270	270 mm	8 mm
A 350008285	285 mm	8 mm
A 350008300	300 mm	8 mm
A 350008315	315 mm	8 mm
A 350008330	330 mm	8 mm
A 350008345	345 mm	8 mm
A 350008360	360 mm	8 mm

A 350008375	375 mm	8 mm
A 350009270	270 mm	9 mm
A 350009285	285 mm	9 mm
A 350009300	300 mm	9 mm
A 350009315	315 mm	9 mm
A 350009330	330 mm	9 mm
A 350009345	345 mm	9 mm
A 350009360	360 mm	9 mm
A 350009375	375 mm	9 mm
A 350010270	270 mm	10 mm
A 350010285	285 mm	10 mm
A 350010300	300 mm	10 mm
A 350010315	315 mm	10 mm
A 350010330	330 mm	10 mm
A 350010345	345 mm	10 mm
A 350010360	360 mm	10 mm
A 350010375	375 mm	10 mm

SOL-T Femur Nail System

- The three cervical screws with an ante-version of 12 deg allow additional stabilisation in the case of ipsilateral femoral neck fractures.
- The optimised anatomical shape (with ante-curvature and lateralisation) allows easy insertion over the tip of the trochanter major.
- Can be used effectively for upper & lower one third femoral fracture applications.



Part No		L Dia (mm)	
Left	Right		
A 345009320	A 346009320	320	9
A 345009340	A 346009340	340	9
A 345009360	A 346009360	360	9
A 345009380	A 346009380	380	9
A 345009400	A 346009400	400	9
A 345009420	A 346009420	420	9
A 345009440	A 346009440	440	9

A 345009460	A 346009460	460	9
A 345009480	A 346009480	480	9
A 345010320	A 346010320	320	10
A 345010340	A 346010340	340	10
A 345010360	A 346010360	360	10
A 345010380	A 346010380	380	10
A 345010400	A 346010400	400	10
A 345010420	A 346010420	420	10
A 345010440	A 346010440	440	10
A 345010460	A 346010460	460	10
A 345010480	A 346010480	480	10
A 345011320	A 346011320	320	11
A 345011340	A 346011340	340	11
A 345011360	A 346011360	360	11
A 345011380	A 346011380	380	11
A 345011400	A 346011400	400	11
A 345011420	A 346011420	420	11
A 345011440	A 346011440	440	11
A 345011460	A 346011460	460	11
A 345011480	A 346011480	480	11



Cervical Screws Dia 6.5 mm

Part No	Length
A 355000060	60 mm
A 355000065	65 mm
A 355000070	70 mm
A 355000075	75 mm
A 355000080	80 mm
A 355000085	85 mm
A 355000090	90 mm
A 355000095	95 mm
A 355000100	100 mm
A 355000105	105 mm
A 355000110	110 mm

SOL-T Instruments Set

Instruments Set for Tibia	A 749000000
Instruments Set for Femur	A 750000000

SOL-T Proximal Femoral Device

- Proximal diameter of the nail 15 mm.
- 6 deg Mediolateral angle of the nail for easy insertion.
- Extended length at Distal end to prevent stress fractures.
- Simple, reliable, versatile & stable fixation.
- User friendly instrumentation.



SOL-T PFD Solid Nail

Part No	Dia
A 360013009	9 mm
A 360013010	10 mm
A 360013011	11 mm



SOL-T PFD Cannulated Nail

Part No	Dia
A 841009242	9 mm
A 841010242	10 mm
A 841011242	11 mm



SOL-T PFD Long Nail System

Part No		L (mm)	Dia
Left	Right		
A 363009340	A 364009340	340	9
A 363009380	A 364009380	380	9
A 363009420	A 364009420	420	9
A 363010340	A 364010340	340	10
A 363010380	A 364010380	380	10
A 363010420	A 364010420	420	10
A 363011340	A 364011340	340	11
A 363011380	A 364011380	380	11
A 363011420	A 364011420	420	11



Hip Bolt 9 mm

Part No	Length
A 361000070	70 mm
A 361000075	75 mm
A 361000080	80 mm
A 361000085	85 mm
A 361000090	90 mm
A 361000095	95 mm
A 361000100	100 mm
A 361000105	105 mm
A 361000110	110 mm
A 361000115	115 mm



Hip Pin 6 mm

Part No	Length
A 362000055	55 mm
A 362000060	60 mm
A 362000065	65 mm
A 362000070	70 mm
A 362000075	75 mm
A 362000080	80 mm
A 362000085	85 mm
A 362000090	90 mm
A 362000095	95 mm
A 362000100	100 mm

SOL-T Proximal Femoral Device Instruments Set A 751000000

Trocantric Femoral Nail



Interlocking Trocantric Nail

Part No	Angle(deg)	Dia
A 412130009	130	9 mm
A 412130010	130	10 mm
A 412130011	130	11 mm
A 412130012	130	12 mm
A 412135009	135	9 mm
A 412135010	135	10 mm

A 412135011	135	11 mm
A 412135012	135	12 mm



Anti Rotational Screw 6.4 mm

Part No	Length
A 413064065	65 mm
A 413064070	70 mm
A 413064075	75 mm
A 413064080	80 mm
A 413064085	85 mm
A 413064090	90 mm
A 413064095	95 mm
A 413064100	100 mm
A 413064105	105 mm
A 413064110	110 mm




Lag Screw 8 mm

Part No	Length
A 414008065	65 mm
A 414008070	70 mm
A 414008075	75 mm
A 414008080	80 mm
A 414008085	85 mm
A 414008090	90 mm
A 414008095	95 mm
A 414008100	100 mm
A 414008105	105 mm

Trocantric Femoral Nail

Instrument set A 921010000

Locking Screws



Locking Screw 3.8 mm

Part No	Length
A 349000020	20 mm
A 349000025	25 mm
A 349000030	30 mm
A 349000035	35 mm
A 349000040	40 mm

A 349000045	45 mm
A 349000050	50 mm
A 349000055	55 mm
A 349000060	60 mm
A 349000065	65 mm
A 349000070	70 mm



Locking Screw 4.5 mm

Part No	Length
A 347000020	20 mm
A 347000025	25 mm
A 347000030	30 mm
A 347000035	35 mm
A 347000040	40 mm
A 347000045	45 mm
A 347000050	50 mm
A 347000055	55 mm
A 347000060	60 mm
A 347000065	65 mm
A 347000070	70 mm
A 347000075	75 mm
A 347000080	80 mm
A 347000085	85 mm
A 347000090	90 mm
A 347000095	95 mm
A 347000100	100 mm



Locking Screw 5 mm

Part No	Length
A 348000025	25 mm
A 348000030	30 mm
A 348000035	35 mm
A 348000040	40 mm
A 348000045	45 mm
A 348000050	50 mm
A 348000055	55 mm
A 348000060	60 mm
A 348000065	65 mm
A 348000070	70 mm
A 348000075	75 mm

A 348000080	80 mm
A 348000085	85 mm
A 348000090	90 mm

3.8mm Screws - To be used for Distal Holes of 8 & 9mm SOL-T Tibia Nails & 9mm SOL-T Femur Nail & of 8mm Ultralock II Tibia Nail.

4.5mm Screws - To be used for Proximal Holes of SOL-T Tibia, Femur Nails when used as conventional Interlock & for Distal Holes of 10 to 12mm Femur & 10mm Tibia Nails. Also used for Proximal Holes of all Ultralock II Tibia Nails and Distal Holes of 9mm Ultralock II Femur and SOL-T, PFD.

5mm Screws - To be used for Proximal Holes of Femur Nail and Proximal and Distal Holes of Supracondylar Nails.

Universal Locking Nail SS



Tibia Nail

Part No	Length	Dia
A 419008270	270 mm	8 mm
A 419008285	285 mm	8 mm
A 419008300	300 mm	8 mm
A 419008315	315 mm	8 mm
A 419008330	330 mm	8 mm
A 419008345	345 mm	8 mm
A 419008360	360 mm	8 mm
A 419008375	375 mm	8 mm
A 419008390	390 mm	8 mm
A 419009270	270 mm	9 mm
A 419009285	285 mm	9 mm
A 419009300	300 mm	9 mm
A 419009315	315 mm	9 mm
A 419009330	330 mm	9 mm
A 419009345	345 mm	9 mm
A 419009360	360 mm	9 mm
A 419009375	375 mm	9 mm
A 419009390	390 mm	9 mm

A 419010270	270 mm	10 mm
A 419010285	285 mm	10 mm
A 419010300	300 mm	10 mm
A 419010315	315 mm	10 mm
A 419010330	330 mm	10 mm
A 419010345	345 mm	10 mm
A 419010360	360 mm	10 mm
A 419010375	375 mm	10 mm
A 419010390	390 mm	10 mm
A 419011270	270 mm	11 mm
A 419011285	285 mm	11 mm
A 419011300	300 mm	11 mm
A 419011315	315 mm	11 mm
A 419011330	330 mm	11 mm
A 419011345	345 mm	11 mm
A 419011360	360 mm	11 mm
A 419011375	375 mm	11 mm
A 419011390	390 mm	11 mm



Femur Nail

Part No	Length	Dia
A 418009320	320 mm	9 mm
A 418009340	340 mm	9 mm
A 418009360	360 mm	9 mm
A 418009380	380 mm	9 mm
A 418009400	400 mm	9 mm
A 418009420	420 mm	9 mm
A 418009440	440 mm	9 mm
A 418009460	460 mm	9 mm
A 418009480	480 mm	9 mm
A 418010320	320 mm	10 mm
A 418010340	340 mm	10 mm
A 418010360	360 mm	10 mm
A 418010380	380 mm	10 mm
A 418010400	400 mm	10 mm
A 418010420	420 mm	10 mm
A 418010440	440 mm	10 mm
A 418010460	460 mm	10 mm
A 418010480	480 mm	10 mm
A 418011320	320 mm	11 mm
A 418011340	340 mm	11 mm

A 418011360	360 mm	11 mm
A 418011380	380 mm	11 mm
A 418011400	400 mm	11 mm
A 418011420	420 mm	11 mm
A 418011440	440 mm	11 mm
A 418011460	460 mm	11 mm
A 418011480	480 mm	11 mm
A 418012320	320 mm	12 mm
A 418012340	340 mm	12 mm
A 418012360	360 mm	12 mm
A 418012380	380 mm	12 mm
A 418012400	400 mm	12 mm
A 418012420	420 mm	12 mm
A 418012440	440 mm	12 mm
A 418012460	460 mm	12 mm
A 418012480	480 mm	12 mm

Universal Locking Nail Instrument Set

For Tibia	AA01110000
For Femur	AA01210000
Tibia with Distal	AA01310000
Femur with Distal	AA01410000

Interlocking Screws



Screw 3.9 mm

Part No	Length
A 421000030	30 mm
A 421000032	32 mm
A 421000034	34 mm
A 421000036	36 mm
A 421000038	38 mm
A 421000040	40 mm
A 421000042	42 mm
A 421000044	44 mm
A 421000046	46 mm
A 421000048	48 mm
A 421000050	50 mm
A 421000052	52 mm
A 421000054	54 mm

A 421000056	56 mm
A 421000058	58 mm
A 421000060	60 mm
A 421000062	62 mm
A 421000064	64 mm
A 421000066	66 mm
A 421000068	68 mm
A 421000070	70 mm



Screw 4.9 mm

Part No	Length
A 420000030	30 mm
A 420000032	32 mm
A 420000034	34 mm
A 420000036	36 mm
A 420000038	38 mm
A 420000040	40 mm
A 420000042	42 mm
A 420000044	44 mm
A 420000046	46 mm
A 420000048	48 mm
A 420000050	50 mm
A 420000052	52 mm
A 420000054	54 mm
A 420000056	56 mm
A 420000058	58 mm
A 420000060	60 mm
A 420000065	65 mm
A 420000070	70 mm
A 420000075	75 mm
A 420000080	80 mm
A 420000085	85 mm
A 420000090	90 mm

3.9mm Screws - To be used for Distal Locking 8mm, 9mm of Universal Tibia & 9mm Universal Femur.

4.9mm Screws - To be used for Proximal locking for Universal Tibia & Femur all diameter nails. Distal locking of TFN all diameter nails.

Interlocking Humerus Nail System



Interlocking Nail For Humerus

Part No.	Dia	Length
A 460006190	6 mm	190 mm
A 460006200	6 mm	200 mm
A 460006210	6 mm	210 mm
A 460006220	6 mm	220 mm
A 460006230	6 mm	230 mm
A 460006240	6 mm	240 mm
A 460006250	6 mm	250 mm
A 460006260	6 mm	260 mm
A 460006270	6 mm	270 mm
A 460006280	6 mm	280 mm
A 460006290	6 mm	290 mm
A 460006300	6 mm	300 mm
A 460006310	6 mm	310 mm
A 460007190	7 mm	190 mm
A 460007200	7 mm	200 mm
A 460007210	7 mm	210 mm
A 460007220	7 mm	220 mm
A 460007230	7 mm	230 mm
A 460007240	7 mm	240 mm
A 460007250	7 mm	250 mm
A 460007260	7 mm	260 mm
A 460007270	7 mm	270 mm
A 460007280	7 mm	280 mm
A 460007290	7 mm	290 mm
A 460007300	7 mm	300 mm
A 460007310	7 mm	310 mm
A 460008190	8 mm	190 mm
A 460008210	8 mm	200 mm
A 460008220	8 mm	220 mm
A 460008230	8 mm	230 mm
A 460008240	8 mm	240 mm

A 460008250	8 mm	250 mm
A 460008260	8 mm	260 mm
A 460008270	8 mm	270 mm
A 460008280	8 mm	280 mm
A 460008290	8 mm	290 mm
A 460008300	8 mm	300 mm
A 460008310	8 mm	310 mm



Interlocking Screw For Humerus Nail Dia 3.3 mm

Part No.	Length
A 459000030	30 mm
A 459000034	34 mm
A 459000036	36 mm
A 459000038	38 mm
A 459000040	40 mm
A 459000042	42 mm
A 459000044	44 mm
A 459000046	46 mm
A 459000048	48 mm
A 459000050	50 mm

Instruments for Interlocking Humerus Nail

Part No. AA00610000



Jig
Part No. A 905001000



Drill Bit Dia 2.7x190 mm
Part No. A 989190027



Drill Sleeve Dia 8 x 2.8 mm
Part No. A 990000000



Tommy Bar
Part No. A 991000000



Awl for Humerus Nail
Part No. A 992000000



Coupling Bolt
Part No. A 993000000



Protection Sleeve 11/8 mm

Part No. A 968000000



Trocar

Part No. A 971000000



Depth Gauge

Part No. A 647001000



Screw Tightener 3.5mm A/F

Part No. A 952001000



Combination Wrench 13 mm

Part No. D 502131000



Extractor Rod

Part No. A 953000000



Extractr Handle

Part No. A 954000000



Extractor Hammer

Part No. A 651000000



Container

Part No. H 934001000

Eight Plate Implants



Eight Plate

Length

12mm

16mm

24mm

32mm

Part No.

A 415000012

A 416000016

A 416000024

A 416000032



Eight Plate Bone Screw

Dia 4.5mm, non cannulated

Length

18mm

20mm

22mm

24mm

26mm

28mm

30mm

32mm

Part No.

A 444000018

A 444000020

A 444000022

A 444000024

A 444000026

A 444000028

A 444000030

A 444000032



Eight Plate Bone Screw

Dia 4.5mm, cannulated

Length

16mm

24mm

32mm

Part No.

A 417000016

A 417000024

A 417000032

Pitkar

Pitkar

Pitkar

Pitkar

Pitkar

Pitkar

Pitkar

Eight Plate Instruments Set

K Wire 1.6 X 150mm

Part No. J 101015016



1.6mm threaded k wire

Part No. AA09201000



Double drill guide

Part No. AA04801000



Cannulated drill bit Dia3.2/1.7

Part No. AA04900000



Direct measuring device

Part No. AA14300000



Self centering drill guide

Part No. AA14401000



A/F 3.5mm solid screw driver

Part No. AA14501000



A/F 3.5mm cannulated screw driver

Part No. AA14601000



Container For Eight Plate Instrument Set

Part No. H988001000

Miscellaneous Implants

K Wires :

Part No.	Length	Size
J100015010	150	1mm
J100015015	150	1.5mm
J100015020	150	2mm
J100015025	150	2.5mm
J100015030	150	3.0mm
J100030010	300	1mm
J100030015	300	1.5mm
J100030020	300	2mm
J100030025	300	2.5mm
J100030030	300	3.0mm

Steinman Pins

Part No.	Size
J102000030	3mm
J102000040	4mm
J102000050	5mm

Endo buttons

Part No. J150000004



SS wire spool

Part No.	Size
BA02100016	16 G
BA02100018	18 G
BA02100020	20 G

